




STERILIZATION INDICATORS



SELF-SEALING POUCHES WITH INTERNAL AND EXTERNAL INDICATORS

The special feature that distinguishes our Self-Sealing Pouches is the fact of having in addition to the usual indicator outside, also an indicator within the welding. The indicator inside is protected by a further welding so that you avoid contact between the material inserted in the roll and the ink as it is required by the regulations.

Benefits of the double parameter of sterilization:

- **Increased safety** during the sterilization process: check done by two indicators, external and internal.
- **Complete guarantee** on the sterilized product: the indicators are in direct contact with the gas that sterilizes the product inserted
- **Saving of money and time:** it is not necessary to insert an indicator strip of sterilization inside the reel/pouch as further check
- **Total and easier check:** regarding the reels, it is possible to check the internal sterilization on all the length of the piece cut
- **Indicators protection:** the internal indicators are protected; they can't enter in contact with the product inserted

A special film with characteristics of resistance and good peelability, sealed to a medical paper with high anti-bacterial protection and good quality to sterilization, are the main components of our flat reels. Raw materials are selected by us, and all the working phases, till to have the finished products, are made under our own direct control in our factory.

FEATURES:

- Medical Paper 60 g/m²
- Laminated Film PET/PP peel blue color
- Indicators STEAM and EO (water based, not toxic)
- Conformity with international standards
- Multi-line weld lines are very durable and equally spaced.
- No release of fibers of paper or lacerations of the plastic material
- Opening direction is marked
- The chevron allows a rapid and correct opening
- In one of the two ends of the envelope has been applied a strip of biadhesive tape to close the pouch

PERSONALIZE YOUR SELF-SEALING POUCHES ON REQUEST:

- Medical Paper 70 g/m² and 80 g/m² available
- Film color available in Green and Transparent
- Add the FORMALDEHYDE indicator
- Print your Private Label on the items
- Personalize also the boxes of the pouches

SELF-SEALING POUCHES

Code	measure in mm
BA60100	60x100
BA75250	75x250
BA90165	90x165
BA90250	90x250
BA14260	140x260
BA19330	190x330
BA19370	190x370
BA19400	190x400
BA30370	300x370
BA30450	300x450
BA30500	300x500